

Note

TRI+ represents the interface between the TRI dental implant system and the exocad CAD/CAM system. The following instructions are intended only for users who are familiar with the exocad system.

Scope

CAD/CAM

- CAD abutments
- Cement retained TRI implant crowns and bridges
- Screw retained TRI implant bars and bridges

Process Description CAD/CAM

Cement Retained Restorations

Cement retained restorations on TRI implants are based on the TRI TV70 titanium bonding bases.


Please use the nt-trading library in which all our titanium bonding bases and scan bodies are implemented.

Screw Retained Restorations

Screw retained restorations on TRI implants are based on the TRI TV40 & TV50 screw retained abutments.


Please import the TRI screw retained implant kit into your exocad software, see download section at the end of the document.

Description:	TV40-01	<input checked="" type="checkbox"/> Connection Geometry
Connection geometry:	Hut-TV40-01offen.stl	<input checked="" type="checkbox"/> Scan Abutment
Scan Abutment:	scanabutment_tri35g_cropped.stl	<input type="checkbox"/> Screw
Screw (optional):		<input checked="" type="checkbox"/> Implant
Implant (optional):	B-TV41B11.stl	<input checked="" type="checkbox"/> Titanium Base
Titanium Base (optional):		
Click Point:	x: 2.92539453 y: 10.1800022 z: 0.00975036	<input type="checkbox"/> Invert connection geometry triangle orientation
		<input type="checkbox"/> Connection geometry has rotation lock ("engaging")
		<input type="checkbox"/> For usage with titanium base





Example: implant kit for TV40-01

The TV40 connection geometry (left picture) is a combination between the TV40 screw retained abutment and the retaining screw RS-PCC. This combination ensures ideal connection geometry and screw channel in the screw retained restoration.




Please consider the following click point for the detection of the scanbody:





Example: implant kit for TV50-34-F



Important Note: the TV50-34-F abutment is asymmetrical due to the angulation. Please use the above specified orientation of the scanbody.

Required Articles from TRI Implants

CADCAM

TV70-07	Titanium Bonding-Base with TRI®-Friction
TV70-20	Titanium Bonding-Base with TRI®-Friction
TV70-Scan	Scanbody for Laboratory and Intraoral Scanner
TV40-01	Screw-Retained Abutment, Straight, 1mm
TV40-02	Screw-Retained Abutment, Straight, 2mm
TV40-04	Screw-Retained Abutment, Straight, 4mm
TV40-06	Screw-Retained Abutment, Straight, 6mm
RS-PCC	Replacement Retaining Screw for PCC
TV50-34-F	Screw Retained Abutment - Angulated 30° - Ø 4.5mm - 4.7mm

Downloads

CADCAM

Please request the following files at digital@tri-implants.com and import into your Exocad system:

- config.xml

For support, please contact: digital@tri-implants.com