

Press Release 10/10/2012

At the 20th EAO convention 2012, TRI Dental Implants Int. AG showcases its innovative TRI®-Friction, SLA surface technology and novel implant packaging.

TRI Dental Implants will showcase its latest innovations as silver sponsor at Europe's biggest implantology congress EAO (European Association of Oral Implantology) in Copenhagen (Denmark) from October 10-13, 2012.

At the core of the showcase is the innovative **TRI®-Friction implant abutment concept**, that can be considered as one of the strongest and best performing in the industry. The **TRI®-Friction** abutment features a one degree taper from the coronal to the apical portion of the outer hex, creating a frictional connection between the implant and abutment. This innovative technology eliminates micro movements and screw-loosening.

Further to this, TRI expands its portfolio with a **SLA surface technology**, which can be considered as gold standard for osseointegration in the industry, and expands the acceptance of wide user groups worldwide.

For the first time, TRI will present its **TRI®-iPack** as a novel solution for implant packaging which will be introduced into the market by the end of this year. The new **TRI®-iPack** guarantees improved handling, better visibility of the implant and storage capabilities.

Company information:

TRI Dental Implants is an innovative Swiss dental implant company. The company follows a valuebased strategy with high-quality innovative products (Swiss made), that consists of an easy-to-use portfolio (50% less components), compatible with leading providers and are positioned in an attractive price bracket (50% under premium providers). TRI Dental Implants was founded in 2010 and is currently distributing to 12 international markets, including high volume markets Germany and Italy. TRI Dental Implants develops and produces its products in Baar, Switzerland.