


Optimized Soft Tissue Management with new Esthetic Healing Collar



TRI® Dental Implants, a young and innovative, swiss high quality provider of dental implants, launches the new, patented gingiva former type „Esthetic Healing Collar“, taking the next step in sustainable improvement of Soft Tissue Management in dental implantology.

The outlined emergence profile of the Esthetic Healing Collar, which is consistently featured in all TRI abutments, has a lasting effect on gingiva growth for optimal esthetic results. As a unique feature our abutment exhibits fixation loops for an exceptionally high surgical attachment of the papillae (see picture).

The Esthetic Healing Collar is available in the emergence profiles 3.5 mm and 5.0 mm for the universal prosthetic platform 3.5 mm with implant diameters of 3.75, 4.1, and 4.7 mm of the TRI-Vent implant system.

Further Information:

Infoline: 00800 3313 3313