

# TRI+




## Note

TRI+ represents the interface between the TRI dental implant system and the CEREC system by Sirona. The following instructions are intended only for users who are familiar with the CEREC system.

## Indications

- CAD abutments
- Cement retained TRI implant crowns and bridges

## Process Description



Required components

Please use the Sirona article **ZTSV3.5 L** for all TRI Vent implants (not compatible with TRI Narrow).

## Downloads

No downloads required.

For support, please contact: [digital@tri-implants.com](mailto:digital@tri-implants.com).